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Technogene Viral Nucleic Acid Extraction Kit

Components of the kit:

- Lysis solution (CN: TG-LS), 1 vial 65 mL
- Si-Nano (CN:SI-NANO), 1vial 10 mL
- Washing solution A (CN: TG-WSA),1 vial 50 mL
- Washing solution B (CN: TG-WSB),1 vial 70mL
- Elution Buffer (CN: TG-DW), 1 vial 10 mL
- High purification column
- Collecton Tube

procedure

1. Add 250 μ L of sample into a microtube

*Note: for Pap Smear sample, transfer as much as possible tissue in to a microtube. after washing with PBS buffer, add 200 μ L of PBS buffer to the pellet and after resuspension , proceed with following steps.

2. Add 600 μ L of Lysis solution to the microtube.
3. Vortex and incubate for 5 mins at 54°C
4. Add 100 μ L Si-Nano solution (Shake Si-nano tube before using)

5. Invert by hand for 2 mins and Transfer lysate to High pure column

*Note: After inverting, vortex the lysate for 3 seconds and immediately transfer to the high pure column, making sure that all nanoparticles are completely transferred.

6. Centrifuge column in 3000 RCF for 3 mins (soft mode).
7. Change the Collection Tube.
8. Add 500 μ L TG Washing solution A to the column and then centrifuge in 5000 RCF for 1 min.
9. Discard transferred Washing solution .
10. Add 700 μ L Washing solution B to the column and centrifuge in 5000 RCF for 1 mins.

11. Discard transferred Washing solution B

- 12.Centrifuge column in upper speed for 3 mins

- 13.Transfer column to new DNase/RNase free microtube .

- 14.Open the cap of the column and incubate in room temperature for 4 mins to drying the column

- 15.Add 50 μ L prewarm(60°C) Elution Buffer and incubate for 3mins

- 16.Centrifuge the column in upper speed for 1min.

Note:to increase yield of Nucleic Acid, you can transfer the product from microtube to column and centrifuge again.

Safety directions: The lysis solution and Washing solution A are tissue hazard material, hence, in the cases of contact with skin, eyes, ETC, wash with water and refer to the hospital for additional medications.

Ratingene.com / Technogene.ir

0913 63 62 862 / 0913 63 62 963

034-32 11 20 45